

P/N: 4233062

Copyright

© 2021, FLIR Systems, Inc.

All rights reserved worldwide. Names and marks appearing herein are either registered trademarks or trademarks of FLIR Systems and/or its subsidiaries. All other trademarks, trade names or company names referenced herein are used for identification only and are the property of their respective owners.

Document identity

Publ. No.: 4233062

Commit: 77402

Language:

Modified: 2021-06-16

Formatted: 2021-06-16

Website

<http://www.flir.com>

Customer support

<http://support.flir.com>

Disclaimer

Specifications subject to change without further notice. Camera models and accessories subject to regional market considerations. License procedures may apply. Products described herein may be subject to US Export Regulations. Please refer to exportquestions@flir.com with any questions.



General description	
A microscope stand with adjustable test article platform compatible with cooled FLIR A and X-series cameras. The microscope stand provides a convenient camera mounting method for testing and IR imaging of objects at a close range.	
Camera mounting adapters are included.	
Physical data	
Weight (without camera)	11 kg (24 lb)
Size (L x W x H) without camera	305 × 305 × 556 mm (12 × 12 × 21.9 in)
Shipping information	
Packaging, type	Cardboard box
Packaging, contents	<ul style="list-style-type: none"> Microscope stand Mounting adapter for cooled A-series cameras Mounting adapter for cooled X-series cameras
Country of origin	USA

Compatible with the following products:

- 29437-261; FLIR A6261 (InGaAs, 180Hz)
- 29440-200; FLIR A6700 (f/2.5, 1-5µm, 60Hz)
- 29440-201; FLIR A6701 (f/2.5, 3-5µm, 60Hz)
- 29440-202; FLIR A6702 (f/4.0, 1-5µm, 60Hz)
- 29440-203; FLIR A6703 (f/4.0, 3-5µm, 60Hz)
- 29440-205; FLIR A6705 (Automation rear panel, 60Hz)
- 29440-250; FLIR A6750 (f/2.5, 1-5µm, 125Hz)
- 29440-251; FLIR A6751 (f/2.5, 3-5µm, 125Hz)
- 29439-251; FLIR A6751 SLS (f/2.5, 7.5-10.5µm, 125Hz)
- 29440-252; FLIR A6752 (f/4.0, 1-5µm, 125Hz)
- 29440-253; FLIR A6753 (f/4.0, 3-5µm, 125Hz)
- 29439-253; FLIR A6753 SLS (f/4.0, 7.5-10.5µm, 125Hz)
- 29440-280; FLIR A6780 (f/2.5, 1-5µm, 125Hz)
- 29440-281; FLIR A6781 (f/2.5, 3-5µm, 125Hz)
- 29440-282; FLIR A6782 (f/4.0, 1-5µm, 125Hz)
- 29440-283; FLIR A6783 (f/4.0, 3-5µm, 125Hz)
- 29439-281; FLIR A6781 SLS (f/2.5, 7.5-11µm, 125Hz)
- 29439-283; FLIR A6783 SLS (f/4.0, 7.5-11µm, 125Hz)



Cooled Science A and X Camera Microscope Stand

P/N: 4233062

© 2021, FLIR Systems, Inc.

#4233062; r. 77402;

- 29443-200; FLIR A8580 (f/2.5, 1.5-5.0 μm)
- 29443-201; FLIR A8581 (f/2.5, 3.0-5.0 μm)
- 29443-202; FLIR A8582 (f/4.0, 1.5-5.0 μm)
- 29443-203; FLIR A8583 (f/4.0, 3.0-5.0 μm)
- 29444-201; FLIR A8581 SLS (f/2.5, 7.5-12.5 μm)
- 29420-200; FLIR X6800sc
- 29420-201; FLIR X6801sc
- 29420-202; FLIR X6802sc
- 29420-203; FLIR X6803sc
- 29267-200; FLIR X6900sc
- 29267-201; FLIR X6901sc
- 29421-201; FLIR X6901sc SLS
- 29267-202; FLIR X6902sc
- 29267-203; FLIR X6903sc
- 29421-203; FLIR X6903sc SLS
- 80601-0101; FLIR X6520sc
- 29422-200; FLIR X8500sc
- 29422-201; FLIR X8501sc
- 29428-201; FLIR X8501sc SLS
- 29422-202; FLIR X8502sc
- 29428-203; FLIR X8503sc SLS
- 29422-203; FLIR X8503sc

ASSEMBLED WEIGHT = 24 LBS

